

BOD WORKSHEET

Example # 1A

	DAY 1	DAY 2	DAY 3	DAY 4	DAY 5
Incubator Temperature °C	20	20.1	20	19.9	20

Location Collected: effluent flume Before Cl₂ ✓ After Cl₂ _____ Sample Seeded: no
 Initial Date/time: 12/14/05 2:00pm Final Date/time: 12/19/05 2:30pm Sample Type: grab Date/Time Sample Collected: 12/13/05 10:00

	INFLUENT (RAW)				EFFLUENT			UP STREAM			DOWN STREAM			BLANKS		
Seeded or <u>unseeded</u>	DUP				DUP			DUP			DUP					
mL used OR % Concentration (P for 100% = 1.0, 75% = 0.75)	15	24	30		75	150	300								300	302
Bottle Number	8	12	27		38	39	40								54	55
Initial DO (D ₁)	8.5	8.4	8.3		8.8	8.8	8.9								8.5	8.5
Final DO (D ₂)	0.5	0.5	0.6		6.9	7.0	7.4								8.4	8.4
Oxygen Demand	8.0	7.9	7.7		1.9	1.8	1.5								0.1	0.1
Seed contribution																
Dilution Factor (1/P)	20	12.5	10		4	2	1									
5-Day BOD (mg/L)	160	98.75	77		7.6	3.6	1.5									

>160 mg/L

<1.5 mg/L

For unseeded samples:

$$\text{BOD}_5, \text{ mg/L} = (D_1 - D_2) \frac{1}{P}$$

None of the dilutions met the 2.0 mg/L minimum depletion. Using 100% sample concentration, multiply 1.5 by the total volume of 300, divided by the highest amount of sample. $<1.5 \times 300/300 = <1.5 \text{ mg/L}$

The detection limit rules at the end of the method deal with how low can you go after you subtract initial – final DO. 3rd bullet at the very end of SM 5210B-20018.b. says detection limit is ~0.1. Because there is a 100% sample dilution in the series, use 1.5 (SM 5210-2001B 7.a.2.) and report as <1.5 mg/L.

Refer to SM 5210B-2001 7.a.3. For diluted influent unseeded, report >160 mg/L (with a greater than as the qualifier due to excessive depletion).

• The lower limit for unseeded samples that require no dilution ($S = 0$; $P = 1.0$) is equal to the detection limit of the DO measurement method (~0.1 mg/L).

3) When all dilutions result in a residual DO < 1.0, select the bottle having the lowest DO concentration (greatest dilution) and report:

$$\text{BOD, mg/L} > \frac{(D_1 - D_2) - (S)V_s}{P}$$

BOD WORKSHEET

Example # 1B

Location Collected: effluent flume Before Cl₂ ✓ After Cl₂ _____ Sample Seeded: no
 Initial Date/time: 12/14/05 2:00pm Final Date/time: 12/19/05 2:30pm Sample Type: grab Date/Time Sample Collected: 12/13/05 10:00

	INFLUENT (RAW)			EFFLUENT			UP STREAM			DOWN STREAM			BLANKS		
Seeded or (unseeded)	DUP			DUP			DUP			DUP					
<u>mL</u> used OR % Concentration (P for 100% = 1.0, 75% = 0.75)	15	24	30		75	150	250							300	302
Bottle Number	8	12	27		38	39	40							54	55
Initial DO (D ₁)	8.5	8.4	8.3		8.8	8.8	8.9							8.5	8.5
Final DO (D ₂)	0.5	0.5	0.6		6.9	7.0	7.4							8.4	8.4
Oxygen Demand	8.0	7.9	7.7		1.9	1.8	1.5							0.1	0.1
Seed contribution															
Dilution Factor (1/P)	20	12.5	10		4	2	1.2								
5-Day BOD (mg/L)	160	98.75	77		7.6	3.6	1.5								

>160 mg/L

<2.4 mg/L

For unseeded samples:

$$\text{BOD}_5, \text{mg/L} = (D_1 - D_2) \frac{1}{P}$$

None of the dilutions met the 2.0 mg/L minimum depletion and there is no 100% sample concentration. Therefore, multiply 2.0 by the total volume of 300, divided by the highest amount of sample (250). $<2.0 \times 300/250 = <2.4 \text{ mg/L}$

The detection limit rules at the end of the method deal with how low can you go after you subtract initial – final DO. 3rd bullet at the very end of SM 5210B-20018.b. says detection limit is ~0.1.

• The lower limit for unseeded samples that require no dilution ($S = 0$; $P = 1.0$) is equal to the detection limit of the DO measurement method (~0.1 mg/L).

Refer to SM 5210B-2001 7.a.3. For diluted influent unseeded, report >160 mg/L (with a greater than as the qualifier due to excessive depletion).

3) When all dilutions result in a residual DO < 1.0, select the bottle having the lowest DO concentration (greatest dilution) and report:

$$\text{BOD, mg/L} > \frac{(D_1 - D_2) - (S)V_s}{P}$$

BOD

Example # 2

Location Collected: effluent flume Before Cl₂ ✓ After Cl₂ _____ Sample Seeded: no
 Initial Date/time: 12/14/05 2:00pm Final Date/time: 12/19/05 2:30pm Sample Type: grab Date/Time Sample Collected: 12/13/05 10:00

	INFLUENT (RAW)			EFFLUENT			UP STREAM			DOWN STREAM			BLANKS		
Seeded or <u>unseeded</u>	DUP			DUP			DUP			DUP					
<u>mL</u> used OR % Concentration (P - 100% = 1.0)	15	24	30		75	150	200							300	302
Bottle Number	8	12	27		38	39	40							54	55
Initial DO (D ₁)	8.5	8.4	8.3		8.8	8.8	8.9							8.5	8.5
Final DO (D ₂)	1.2	1.1	1.0		6.9	7.0	7.4							8.4	8.4
Oxygen Demand	7.3	7.3	7.3		1.9	1.8	1.5							0.1	0.1
Seed contribution															
Dilution Factor (1/P)	20	12.5	10		4	2	1.5								
5-Day BOD (mg/L)	146	91.2	73		7.6	3.6	2.25								

>146 mg/L

>7.6 mg/L

For unseeded samples:

$$BOD_5, \text{ mg/L} = (D_1 - D_2) \frac{1}{P}$$

None of the dilutions met the 2.0 mg/L minimum depletion. There's evidence of toxicity in sample. Report as (greater than) > 7.6 mg/L

All the dilutions met the 1.0 mg/L final DO. However, there appears to be some toxicity in the influent. SM 5210B 7.b. (see below) So, use largest dilution and report as (greater than) >146.

Samples showing large differences between the computed BOD for different dilutions, for example, greater than 30%, may indicate the presence of a toxic substance or analytical problems. When the effect becomes repetitive, investigate to identify the cause. Identify results in the test reports when any of

Example # 3

Location Collected: effluent flume Before Cl₂ ✓ After Cl₂ _____ Sample Seeded: no
 Initial Date/time: 12/14/05 2:00pm Final Date/time: 12/19/05 2:30pm Sample Type: grab Date/Time Sample Collected: 12/13/05 10:00

	INFLUENT (RAW)				EFFLUENT			UP STREAM			DOWN STREAM			BLANKS		
				DUP				DUP			DUP			DUP		
mL used OR % Concentration (P)	15	24	30		75	150	300								300	302
Bottle Number	8	12	27		38	39	40								54	55
Initial DO (D ₁)	8.5	8.4	8.3		8.8	8.8	8.9								8.5	8.5
Final DO (D ₂)	4.5	2.4	0.8 X		0.5	0.5	0.5								8.4	8.4
Oxygen Demand	4.0	6.0	7.5		8.3	8.3	8.4								0.1	0.1
Seed contribution																
Dilution Factor (1/P)	20	12.5	10		4	2	1									
5-Day BOD (mg/L)	80	75	X		31											

Avg. 78 mg/L > 31 mg/L

For unseeded samples: BOD₅, mg/L = $(D_1 - D_2) \frac{1}{P}$

Refer to 5210B-2001 7.a.3.

3) When all dilutions result in a residual DO < 1.0, select the bottle having the lowest DO concentration (greatest dilution) and report:

$$\text{BOD, mg/L} > \frac{(D_1 - D_2) - (S)V_s}{P}$$

Where S = 0 for unseeded samples

The effluent is a bit tricky, but since final DO < 1, set final DO at 1 and use greatest dilution (lowest DO concentration).

8.8 - 1 = 7.8 x 4 = greater than (>) 31 (rounded to whole value). The greater than (>) sign is used as the qualifier due to excessive depletion.

where:

- D₁ = DO of diluted sample immediately after preparation, mg/L,
- D₂ = DO of diluted sample after 5 d incubation at 20°C, mg/L,
- S = oxygen uptake of seed, Δ DO/mL seed suspension added per bottle (¶ 6d) (S = 0 if samples are not seeded),
- V_s = volume of seed in the respective test bottle, mL, and
- P = decimal volumetric fraction of sample used; 1/P = dilution factor.

Example # 4

	DAY 1	DAY 2	DAY 3	DAY 4	DAY 5
Incubator Temperature °C	20	20	20	20	20

Location Collected: inf/eff flows Initial Date: 3-10-05 13:16 After Cl₂: ✓ Sample Seeded: ✓
 Sample Temp: 4°C Final Date: 3-15-05 14:28 Before Cl₂: _____ Time Sample Collected: 07:30

	INFLUENT (RAW)				EFFLUENT				UP STREAM			DOWN STREAM			BLANKS		
	2 mL SEED			DUP	2 mL SEED			DUP	DUP			DUP			300ml	300ml	
mL. used OR % Concentration (P)	6	9	12		150	200	300	300								100%	100%
Bottle Number	326	337	292		307	311	312	313								8.14	8.16
Initial DO (D ₁)	8.09	8.02	7.93		8.89	8.91	8.98	8.97								7.97	8.00
(-) Final DO (D ₂)	4.39	2.73	0.72	X	7.16	7.09	7.04	7.02								0.17	0.16
Oxygen Demand (D ₁ - D ₂)	3.70	5.29	7.21		1.73	1.82	1.94	1.95									
(-) Depletion due to seed (B ₁ - B ₂)f	1.66	1.66	1.66		1.66	1.66	1.66	1.66									
Oxygen Demand after seed correction	2.04	3.63	5.55		0.07	0.06	0.28	0.29									
Dilution Factor (1/P)	50	33.3	25		2.0	1.5	1.0	1.0									
5-Day BOD (mg/L)	102	120.9	X														

See next page for
seed calculations

Avg. 111 mg/L Avg. 0.29 mg/L

For unseeded samples: $BOD_5, \text{ mg/L} = \frac{(D_1 - D_2)}{P}$ For seeded samples: $BOD_5, \text{ mg/L} = \frac{(D_1 - D_2) - (B_1 - B_2) f}{P}$

For the 4th case, BOD, only 1mL of seed should have been used; 2mL is too high. But, otherwise it appears the green influent is calculated correctly. Average the 2 values meeting depletion criteria before rounding to whole number.

The effluent is tricky, but because you have a 100% sample (5201B 7.a.2.), you can use actual depletion even though <2, so sample would be 0.29.

2) If DO depletion is less than 2.0 mg/L and sample concentration is 100% (no dilution except for seed, nutrient, mineral, and buffer solutions), actual seed-corrected, DO depletion may be reported as the BOD even if it is less than 2.0 mg/L.

	Seed Control			
mL. used OR % Concentration (P)	5	7	9	Avg
Initial DO (B ₁)	8.15	8.17	8.16	
Final DO (B ₂)	3.57	2.42	1.32	
Oxygen Demand	4.58	5.75	6.84	
(B ₁ -B ₂)f	1.83	1.64	1.52	1.66
Oxygen Demand per mL of seed material	4.58/5=0.92	5.75/7=0.82	6.84/9=0.76	Avg 0.83

$$f \text{ Value: } f = \frac{(\text{volume of seed in sample dilutions})}{(\text{volume of seed in seed control})}$$

✓ The DO uptake attributable to seed added to each bottle generally should be between 0.6 and 1.0 mg/L. However, the seed added should be adjusted from this range to that required to provide acceptable GGA results.

Glucose-Glutamic Acid (GGA) Test

5210B-2001

Intended use for evaluation of dilution water quality, seed effectiveness, and analytical technique

Use 6 mL of BOD Standard and 2 mL of seed (typically 2 mL)

	GGA std			
mL GGA used	6	6		
mL seed	2	2		
Bottle Number	296	297		
Initial DO (D ₁)	8.21	8.22		
(-) Final DO (D ₂)	2.15	2.32		
Oxygen Demand (depletion) (D ₁ - D ₂)	6.06	5.90		
(-) Depletion due to seed (B ₁ - B ₂)f	1.66	1.66		
Net depletion due to GGA (D ₁ - D ₂) - (B ₁ - B ₂)f	4.40	4.24		
Dilution Factor (300/6)	50	50		
5-Day BOD (mg/L)	220	212		

Avg. 216 mg/L

GGA std BOD values should lie within the range of 198 ± 30.5 mg/L. (167.5-228.5 mg/L)

1. Weak seed usually causes LOW GGA results.
2. Some sewage seeds are relatively inactive and yield LOW GGA results.
3. Soap contamination typically yields HIGH GGA results.
4. Distilled water contaminated with copper yields LOW GGA results.
5. Too much seed usually causes HIGH GGA results
6. BOD bottles that are not properly rinsed is indicated by HIGH GGA results and would also result in high blank depletion