## Submitting Specimens for Measles Immunity Testing For Exposed Individuals

## **Collection**

- Whole blood collected in a red top tube
- Ideally, 2 mLs should be collected, but the test can be run on less
- Each tube needs to have two patient identifiers on it

## **Storage and Submission**

- Specimen should be stored at 2-8°C
- Ship on cold packs within 48 hours
  - If specimen cannot be shipped within 48 hours:
    - The specimen will need to be spun down, aliquoted, frozen, and shipped on dry ice
- A PH-4182 (Clinical Requisition Form) will need to be completed for each specimen
  - Under Test Requested check Other Miscellaneous and write in Measles IgG testing (see below)
- There is no option for Measles IgG in PTBMIS; an electronic order should **not** be submitted
- Please send a line list of specimens being sent to the SPHL to Amanda Hartley at VPD.Imm@tn.gov.
- Send by courier to the Public Health Lab

*Test Requested		·
Culture	Parasitology	<u>Molecular</u>
□ Actinomycete (Aerobic)	□ Blood Parasite	□ Bordetella (Pertussis) PCR
□ Aerobe	□ Ova & Parasite	□ C.trachomatis/N.gonorrhoeae (GenProbe)
□ Anaerobe	□ Cryptosporidium	□ GI Panel (Biofire)
□ Enteric	Serology	□ ESBL
□ Gonorrhea	□ Arbovirus Panel	□ Norovirus PCR
□ Herpes Simplex Virus	□ HBV Screen	□ Plasmodium PCR
□ Legionella	☐ HCV Screen	ARLN
□ Mycobacteria Smear & Culture	☐ HIV Screen	C. auris Colonization
□ Mycobacteria Reference Isolate	□ Measles/Rubella IgM	□ Candida species Confirmation
□ Mycology	□ Syphilis RPR	□ CRE/CRPA/CRAB Colonization
□ Viral: Virus Suspected	□ VDRL	□ CRE/CRPA/CRAB Confirmation
Other Miscellaneous (Please specify)	Measles IgG testing	