

Division of Laboratory Services

630 Hart Lane Nashville, TN 37216 615-262-6300

https://www.tn.gov/newbornscreeninglab

Newborn Screening Disorder:

Severe Combined Immunodeficiency (SCID)

Alternate Name(s)	 Severe Combined Immunodeficiency, X-Linked,T Cell-Negative, B Cell-Positive, NK Cell-Negative SCIDX XSCID SCID, X-LINKED Immunodeficiency
Analyte(s) Tested	 T Cell Receptor Excision Circles (TREC) Ribonuclease P30 (RPP30)
Methodology	Real-Time Quantitative Polymerase Chain Reaction (RT-qPCR)
TDH Requisition Form	 PH-1582 Form Requests: Contact state lab by email or fax. Include facility name, address, phone number and contact person on your request. Fax: (615) 262-6455 Email: DCLAB.supply@tn.gov
Acceptable Specimen	Dried blood spots on filter paper collected from infants less than 6 months of age
Collection Information	 Optimal specimen: Collect at 24 hours + 1 minute of life Acceptable specimen: Collect 24-48 hours of life If transfused: Recollect 4 days post transfusion
Shipping Information	Health Departments and Birthing Hospitals Private Clinics and Midwives
Screening Results	Tennessee Newborn Screening's Secure Remote Viewer (SRV) Healthcare providers must be registered to view and print patient result reports. To Register: Complete and submit SRV Access Form (PH-3909). NBS List of Screened Disorders and Mailer Explanations
Laboratory Location Performing Testing	Nashville, TN